

EMP3-epithelial membrane protein 3 Gene

Reference: PTXBC009718

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	YMP; epithelial membrane protein 3; EMP-3; HNMP-1; hematopoietic neural membrane protein 1		

For research use only.

GENE DESCRIPTION	epithelial membrane protein 3
GENE ID	2014
ORIGIN SPECIES	Human
ACCESSIONS	BC009718
VECTOR	pDONR223
NCBI SYMBOL	EMP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcactcctctgtgtgtgtcagccctcacatcctcattcttatactgcttttctggccactttggacaagtcctgtgtgactctccctgggaaagatccctgaatctctgtgtacactgcacgtggaacaacacacaaacatgggcctgcagtaatgtcagcgagaatggctgg
ctgaaggcgtgtcaggtcctcatggtgtctcctcattctctgtgtctctcctcattctgtgtctctcctcctgtgtcatgttccagctctacacatgcgacgaggaggtctctctatgccaccggcctctgccagctttgcaccagcgtggcgtgttttactgtgcgcttgatctatgccatccacgcgaggaga
tcctggagaagcaccgcgagggggcagcttcggatactgtctgccttgccctgggtggccttccccctgcctgtgtcagcggcatcatctacatccactacggaagcgggagtg