

MST150-MSTP150 Gene

Reference: PTXBC009719

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MST150; C5orf62; small integral membrane protein 3; NGF-induced differentiation clone 67 protein; small membrane protein NID67		

For research use only.

GENE DESCRIPTION	MSTP150
GENE ID	85027
ORIGIN SPECIES	Human
ACCESSIONS	BC009719
VECTOR	pDONR223
NCBI SYMBOL	MST150
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcagtcagccaagtcccatggaagtcgtgctccaagcacatcctgatatctgggtattgtcctcatcatcctggccaccattgtcatcatgacctcgtttgtgtgtgccagccactgcagtaatcatctatgcgatgcggactcatccgatccttagtggggctgtttga