

TP53RK-TP53 regulating kinase Gene

Reference: PTXBC009727

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	atypical serine/threonine protein kinase TP53RK; EKC/KEOPS complex subunit TP53RK; BUD32; C20orf64; Nori-2; Nori-2p; dJ101A2; TP53-regulating kinase; p53-related protein kinase; TP53 regulating kinase		

For research use only.

GENE DESCRIPTION	TP53 regulating kinase
GENE ID	112858
ORIGIN SPECIES	Human
ACCESSIONS	BC009727
VECTOR	pDONR223
NCBI SYMBOL	TP53RK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccagagctactacgccggccgatggcaggagcccccggagctgaggtctggccgcagccgggagcggagcagccgtcttgagcggcctggagctgggaagcagggtgccgaggcgcgctgttcggtgccgcttcaggcgccgcggcggtgatcaagcaccgcttcccaagggtaccgacccggcgctggagggcgggctggcagacggcgagcgggtgagggagccggcgctcccgctgctgcggcgctgggaatctgccccagttgtctttttgtggactatgcttccaactgcttatatggaagaattgaaggctcagtgactgttcgagattatattcagtcactatggagactgaaaaaacccccagggtcttccaacttagccaagacaattgggcaggttttggctcgaatgcacgatgaagacctcattcatggtgatctcaccacctccaacatgctctgaaacccccctggaaacagctgaacattgtgctcatagactttggcgtgagtttcatcttcagcactccagagataaaggagtagacctctatgctcctggagaaggccttctcagtagccatcccaacactgaaactgtgttgaagcctttctgaagagctactccacctcccaaaaggccaggccagtgctaaaaaaattagatgaagtgcgcctgagagggaagaaagaggtccatggttggtag