



Reference: PTXBC009748

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBDR; DRD1B; DRD1L2; D(1B) dopamine receptor; D1beta dopamine receptor; d(5) dopamine receptor; dopamine D5 receptor; dopamine receptor D1B; dopamine receptor D5		

For research use only.

GENE DESCRIPTION	dopamine receptor D5
GENE ID	1816
ORIGIN SPECIES	Human
ACCESSIONS	BC009748
VECTOR	pDONR223
NCBI SYMBOL	DRD5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]