



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TRIM41-tripartite motif-containing 41 Gene

Reference: PTXBC009762

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase TRIM41; RINCK; RING finger-interacting protein with C kinase; RING-finger protein that interacts with C kinase; tripartite motif containing 41		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 41
GENE ID	90933
ORIGIN SPECIES	Human
ACCESSIONS	BC009762
VECTOR	pDONR223
NCBI SYMBOL	TRIM41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgcggatgctctgcagcctgcgagagtggaactgacgttggaacctgcacacggctaccacggccctgactgctgccctgaccgcgggggggtccgcctgagagacggcgacgagagtgctgacacatccaagcctctcggccgactgctgctactgggggcccag
ggcttcgcctcgcggccgacatactggagagtagaggtggggggggcgccggggctggggctgggtgctgccctgtaatacaacccatcataaggaanaaggtggggccctgggggttcctcctggggcagcggggatgccagctcctcgcgcacacatcgcggccggcgctcc
actcggccacagccccctgcctcagcgggaagtgtgtgctggtggggcacaacacgcaaacgctatcagccggacgagcctccacagaacagacgctcgtgagccccaagtgaanaaaccaacgctgttggctgtgtaacctgacatgaagatcgccgcctggcctctacacacgagaga