

AGK-acylglycerol kinase Gene

Reference: PTXBC009775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CATC5; CTRCT38; MTDPS10; MULK; acylglycerol kinase, mitochondrial; hAGK; hsMuLK; multi-substrate lipid kinase; multiple substrate lipid kinase		

For research use only.

GENE DESCRIPTION	acylglycerol kinase
GENE ID	55750
ORIGIN SPECIES	Human
ACCESSIONS	BC009775
VECTOR	pDONR223
NCBI SYMBOL	AGK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggtgttctttaaagcgttcgaaatcactggaagaaaactacagctgggctctgcctgctgacctggggaggccattggctctatggaaaacactgtgataacctcctaaggagagcagcctgtcaagaagctcagcactatcaggatgaatcacgctgggagccaactctgagtag
aactccaggaggttaa