

ODF2L-outer dense fiber of sperm tails 2-like Gene

Reference: PTXBC009779

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ977L11.1; outer dense fiber protein 2-like; outer dense fiber of sperm tails 2 like		

For research use only.

GENE DESCRIPTION	outer dense fiber of sperm tails 2-like
GENE ID	57489
ORIGIN SPECIES	Human
ACCESSIONS	BC009779
VECTOR	pDONR223
NCBI SYMBOL	ODF2L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgctagaaaatctaactgataatgaaagtgaataactaacttaagaagaaggattttgaaaaggaggcccatccaagaacttctgtttgagagtgaaaagagcttagaaaccaagatagccaagtggaaacctgcaatcaagaatgaataagaatgaggctatagtgatgaa
agaagcaagtaggcaaaaaactgtagctttaaaaaaggcatctaa