

MRPS14-mitochondrial ribosomal protein S14 Gene

Reference: PTXBC009788

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DJ262D12.2; HSMRPS14; MRP-S14; S14mt; 28S ribosomal protein S14, mitochondrial; mitochondrial 28S ribosomal protein S14; mitochondrial ribosomal protein S14		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S14
GENE ID	63931
ORIGIN SPECIES	Human
ACCESSIONS	BC009788
VECTOR	pDONR223
NCBI SYMBOL	MRPS14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcctcatgctgggctgctgctgggacgttcaagcagatggttcctcatcagcttcaggccaagttcgaagtcactatgtagactggagaatgtggcgcgatgtgaagagacgaaatggcctatgaatacgagatgagggctacgtattaattcactcaggaagaatac
cattttgccaaaaattctcaggatgtggctgatgaagaaatgctgccctccccgggtagctgtcctgttagaatcagaatacgggtgtgtatgacgtccccgtccgctggtgtgaagcggcgcctggaggcttagctgtatagctctccgtcacttagctgacctgggcaacttctgggac
cagcgagcgacatggtaa