

DAD1-defender against cell death 1 Gene

Reference: PTXBC009798

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	oligosaccharyl transferase subunit DAD1; dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1; OST2; DAD-1; oligosaccharyltransferase 2 homolog; oligosaccharyltransferase subunit 2 (non-catalytic); defender against cell death 1		

For research use only.

GENE DESCRIPTION	defender against cell death 1
GENE ID	1603
ORIGIN SPECIES	Human
ACCESSIONS	BC009798
VECTOR	pDONR223
NCBI SYMBOL	DAD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgctcgttagtctgtcatttcgcgggttctagaagagtacttgagctccattccgcagcgtctgaagttgctggacgcgtacctgtgtatatactgctgaccggggcgctgcagttcgggttactgtctcctcgtggggaccttcccttcaactctttctctcgggcttcattctgtg
tggggagtttcatcctagcgggttgcctgagaatacagatcaaccacagaacaagcggaattccaaggcatctccccagagcagcctttgctgatttctctttgccagcaccatctgcacctgtgtgtcatgaactttgttgctga