



Reference: PTXBC009808

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NBR1, autophagy cargo receptor; 1A1-3B; IAI3B; M17S2; MIG19; next to BRCA1 gene 1 protein; B-box protein; cell migration-inducing gene 19 protein; membrane component, chromosome 17, surface marker 2 (ovarian carcinoma antigen CA125); migration-inducing pr		

For research use only.

GENE DESCRIPTION	neighbor of BRCA1 gene 1
GENE ID	4077
ORIGIN SPECIES	Human
ACCESSIONS	BC009808
VECTOR	pDONR223
NCBI SYMBOL	NBR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr