

BIN3-bridging integrator 3 Gene

Reference: PTXBC009824

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bridging integrator 3		
For research use only.			

GENE DESCRIPTION	bridging integrator 3
GENE ID	55909
ORIGIN SPECIES	Human
ACCESSIONS	BC009824
VECTOR	pDONR223
NCBI SYMBOL	BIN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggattcctttaaagattggcagcccaagaacagattgtgcccaaacagtgaggagagactttgaaagggagtgatgaaacttcagcagctggaagagcagaccggaggctgcagaagacatgaagaagagcaccgacgacagctggccatgtcaaatctgccg
tgaagatccttggacttactctccatccctctgtgagcaagaccaggaccttgaacatggtagccctggacacggccatgaagcggatggatgcttcaatcaggaaggtgaaccagatccagaagactgtgatcgaagccctaaaaaagttcggcagtgcttcccagc
ctcaacatggctgtgaagaggcgggaacagccttgcaggactacaggagctgcaggccaaggtggagaagtatgaggaagagagaagacggggccagtgtctggccaagctccaccagcagcagagaggagctgcggcctgtcgggaggacttgaagccaagaacagc
agctgctggaggagatgccgcgttctacggcagccgcctcgaacttccagccagctttagtccctatccgagctcaggttgtgtactactcggaaatgcacaagatcttggagacctgtccatcagcttgaccagccaggccactccgatgagcagcgggagcgggagaacg
aggccaactcagtgagctccgggcctctcattgtggccgatgactga