

RDH13-retinol dehydrogenase 13 (all-trans/9-cis) Gene

Reference: PTXBC009881

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SDR7C3; retinol dehydrogenase 13; retinol dehydrogenase 13 (all-trans and 9-cis); retinol dehydrogenase 13 (all-trans/9-cis); short chain dehydrogenase/reductase family 7C member 3		

For research use only.

GENE DESCRIPTION	retinol dehydrogenase 13 (all-trans/9-cis)
GENE ID	112724
ORIGIN SPECIES	Human
ACCESSIONS	BC009881
VECTOR	pDONR223
NCBI SYMBOL	RDH13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagtgtgaggcgcagcaaaaggacatccgcggggagaccctcaatcaccatgtcaacgcccggcactggactggctccctcaagtctatccgagagtttcagcaaaagatcattgaagaggaggagcagtgacattctaatacaaacgcgggtgtgatgcgggtgcc
cccactggaccaccgaggacggcttcgagatcagtttggcgttaaccacctgggtcactttctcttgacaaacttgctgctggacaagctgaaagcctcagccccttcgcggatcatcaacctctctccctggccatgttctgggcacatagactttgacgacttgaactggcagacga
ggaagtataacaccaaagccgctactgccagagcaagctcgcctctctcaccaggagctgagccggcgctgcaaggctctgtgtgactgtcaacgccctgcacccggcgtggccagagacagctgggcagacacacgggcatccatggctccacttctcagca
ccacactcgggcccatctctgtgctgtgtcaagagccccgagctggccgcccagcccagcacatacctggcgtggcggaggaaactggcggatgttccgaaagtacttcgatggactcaacagaaggccccggccccgaggtgaggtgaggtggcccgaggctt
tggcgtgaaagtgcggcctgtgggttagaggctccctctgtgaggggagcagccctccccagataa