



Reference: PTXBC009911

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GID2; GID2B; protein RMD5 homolog B; GID complex subunit 2 homolog B; required for meiotic nuclear division 5 homolog B		

For research use only.

GENE DESCRIPTION	required for meiotic nuclear division 5 homolog B (<i>S. cerevisiae</i>)
GENE ID	64777
ORIGIN SPECIES	Human
ACCESSIONS	BC009911
VECTOR	pDONR223
NCBI SYMBOL	RMND5B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]