



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FRAG1-FGF receptor activating protein 1 Gene

Reference: PTXBC009930

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

FRAG1; CWH43-N; HPMRS3; MRT17; MRT21; post-GPI attachment to proteins factor 2; FGF receptor activating protein 1; cell wall biogenesis 43 N-terminal homolog; mental retardation, non-syndromic, autosomal recessive, 21; post-GPI attachment to proteins 2

For research use only.

GENE DESCRIPTION	FGF receptor activating protein 1
GENE ID	27315
ORIGIN SPECIES	Human
ACCESSIONS	BC009930
VECTOR	pDONR223
NCBI SYMBOL	FRAG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtaccaggctccaccaccactggatggcggtgagacctcgtacggctcgccttaccatggtggccctgagcaggtctcgttcacattctgccttctctctctgcatcctctggtcctctcttcacactcaaggagacaacgccacacactgtgggggtcccaattaccctccctc
gttgagctcagccatcggcggggaggtgcccagcgctacgtgtggcggtttctgcatcggcgtcactcggcgctcgtcttgggtggcctcgtcctactggaaccatacctcagctgcacctccccgtgtctctgctatcgccgcctctccgcctcaacttgcggctcaatgtctgtga
gaacctcgcgtgtcctagtctcactatgtctctctccaggagcttaccatccacagaaaagcttcaatgtgtgttcatctgcctcatccctcgggacatgctcctcactctcattctcggcggttgaccaagaagcacacagatgcgaagtctctacagctggaacaagcgctcttcatcatc
aatctatctctctctctcggcgtcgtctctatttgcgcacaacatgtattgtgagcgctggagctgacacattcttgcacatcctggagtagactgtgtcttaaccaacatggcgttccacatgacggcctgtggggacttcgggaacaaggagctgctcataacctctcagctcaggaaaa
agcattctca

More info: Please contact us for any question: contact@proteogenix.fr