



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KCTD17-potassium channel tetramerisation domain containing 17 Gene

Reference: PTXBC009961

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

BTB/POZ domain-containing protein KCTD17; potassium channel tetramerization domain containing 17

For research use only.

GENE DESCRIPTION	potassium channel tetramerisation domain containing 17
GENE ID	79734
ORIGIN SPECIES	Human
ACCESSIONS	BC009961
VECTOR	pDONR223
NCBI SYMBOL	KCTD17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr