



Reference: PTXBC009967

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GKLP; HT019; NKTLL; NTKL; P105; SCAR21; TAPK; TEIF; TRAP; N-terminal kinase-like protein; SCY1-like protein 1; SCY1-like, kinase-like 1; coated vesicle-associated kinase of 90 kDa; likely ortholog of mouse N-terminal kinase-like protein; telomerase regulat		

For research use only.

GENE DESCRIPTION	SCY1-like 1 (<i>S. cerevisiae</i>)
GENE ID	57410
ORIGIN SPECIES	Human
ACCESSIONS	BC009967
VECTOR	pDONR223
NCBI SYMBOL	SCYL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr