



Reference: PTXBC009999

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MUS81 structure-specific endonuclease subunit; MUS81 endonuclease homolog; crossover junction endonuclease MUS81; SLX3; SLX3 structure-specific endonuclease subunit homolog		

For research use only.

GENE DESCRIPTION	MUS81 endonuclease homolog (S. cerevisiae)
GENE ID	80198
ORIGIN SPECIES	Human
ACCESSIONS	BC009999
VECTOR	pDONR223
NCBI SYMBOL	MUS81
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr