



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

C8orf33-chromosome 8 open reading frame 33 Gene

Reference: PTXBC010001

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0488 protein C8orf33; chromosome 8 open reading frame 33		
For research use only.			

GENE DESCRIPTION	chromosome 8 open reading frame 33
GENE ID	65265
ORIGIN SPECIES	Human
ACCESSIONS	BC010001
VECTOR	pDONR223
NCBI SYMBOL	C8orf33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]