

MRAP2-melanocortin 2 receptor accessory protein 2 Gene

Reference: PTXBC010003

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf117; bA51G5.2; melanocortin-2 receptor accessory protein 2; MC2R accessory protein 2; melanocortin 2 receptor accessory protein 2		

For research use only.

GENE DESCRIPTION	melanocortin 2 receptor accessory protein 2
GENE ID	112609
ORIGIN SPECIES	Human
ACCESSIONS	BC010003
VECTOR	pDONR223
NCBI SYMBOL	MRAP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgccagaggttaatttctacagaacctccagcaatcgccatcattctgattacacctgggaatgaatattatgagattggaccagttctcttgaaggactgaaggtcataaatattccattgtgattggatttgggttggctctgcagcttcgtgattttatgtttttgtgctgac
cttgcctaccaagacagggagccccacccaagacaatgcagagtcctcagagaagagattcagaatgaacagcttctgtcagacttgggaagacctctggagccagataaagtattttctgccaaaggcaacgagagagtcaggctctctcttctactgctacatcaatgaggtggaacgctt
ggacagagccaagcttgcaccagaccagcccttgacagtgacgtccaactccaggaagccatcagaagcagtgggcagccagaggagctgaacaggtcatgaagttgacatcccaacttctgaacacagaccagaactactttggggagagatgactctctgattctga
accacattattctctgaaactaagccacttccagacctcacaaagacctgattga