



Reference: PTXBC010005

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf67; PPP1R121; phosphorylated CTD-interacting factor 1; PDX-1 C terminus-interacting factor 1; protein phosphatase 1, regulatory subunit 121; PDX1 C-terminal inhibiting factor 1		

For research use only.

GENE DESCRIPTION	PDX1 C-terminal inhibiting factor 1
GENE ID	63935
ORIGIN SPECIES	Human
ACCESSIONS	BC010005
VECTOR	pDONR223
NCBI SYMBOL	PCIF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr