

WBP1-WW domain binding protein 1 Gene

Reference: PTXBC010012

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WBP-1; WW domain-binding protein 1; WW domain binding protein 1		
For research use only.			

GENE DESCRIPTION	WW domain binding protein 1
GENE ID	23559
ORIGIN SPECIES	Human
ACCESSIONS	BC010012
VECTOR	pDONR223
NCBI SYMBOL	WBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcggccagcagcgggaacggcagcagggagcctgggggcaacttcggcgccgcaacagcagctcgagagctgtgccaggagtgacaaccagccctacctgtgagagtggtcactgctcggggagactggctgctgcacctactatgagctctggtggtt
ctggctgctctggactgtctcctcctctttagctgtgtgctgctccgccaccgacgagctaaactcaggtgcacaacacagcagcggcagcgtgaaatcaactgttggcctatcatggggcatgccatgggctggtccttccctaccgggttactgcttgacctgcttctcagcac
ctcaagccccagcctacgaggtatgtgttcaccgcccagcgcacacaccccccttatactgtggcccagcggcccttgactgcttccagtgaacaacctgctgttccctcctaccagctgccctgccactttgaaggacaatgtggaagggtttctccaccagagtc
ccccctcatcagaggggtgagccggggcaggggtgacccctgcctccacacccccctctgcccgtatcggcgtttaaactggcgactccggtattgagctctgcccctgtcctgctcccggtgagggtgagccagtcaggaggtgagggttagtgccacctgccagatctggagg
actactccccgtgtgcactacccccagagctgtaccgcagatcttccatgggctgtcttcagtgaggggacatccataa