

DRAP1-DR1-associated protein 1 (negative cofactor 2 alpha) Gene

Reference: PTXBC010025

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NC2-alpha; dr1-associated corepressor; negative co-factor 2-alpha; negative cofactor 2 alpha; DR1 associated protein 1		

For research use only.

GENE DESCRIPTION	DR1-associated protein 1 (negative cofactor 2 alpha)
GENE ID	10589
ORIGIN SPECIES	Human
ACCESSIONS	BC010025
VECTOR	pDONR223
NCBI SYMBOL	DRAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgagcaagaagaagtaaacgcgcgttcccgccggcgcgatcaagaagatcatgcagacggacgaagattgggaaggtggcggcgccgtgcctgtcatcatctccggcgctcgagctctcctagagtcgctgttgaagaagcctgccaggtgacccagtc
gcggaacgcgaagacatgaccacatccacctgaagcagtgcatcgagctggagcagcagtttgaactcttgaaggacctggtgcatctgttcccgacatgcagggggacggggaagacaaccacatggatggggaagggcgcccgagggccggaagccagcgagcgg
cggccggaagaacggtgggatgggaacaaaagcaaggacaagaagctgtccgggacagactcggagcaggaggatgaatctgaggacacagatactgatggggaagaggagacatcacacccccacccagcagccacccctctgccactttcagagcccccgacac
ccttctgcccttcgctctactctgcctttgccccagcgcggccggccctcagcacctgatgaagaggacgaagaagattacgactcctag