



Reference: PTXBC010042

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IMP4 homolog, U3 small nucleolar ribonucleoprotein; U3 snoRNP protein IMP4; IMP4, U3 small nucleolar ribonucleoprotein, homolog; IMP4, U3 small nucleolar ribonucleoprotein; U3 small nucleolar ribonucleoprotein protein IMP4; BXDC4; U3 snoRNP protein 4 homo		

For research use only.

GENE DESCRIPTION	IMP4, U3 small nucleolar ribonucleoprotein, homolog (yeast)
GENE ID	92856
ORIGIN SPECIES	Human
ACCESSIONS	BC010042
VECTOR	pDONR223
NCBI SYMBOL	IMP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]