



Reference: PTXBC010047

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 SUMO-protein ligase PIAS4; PIAS-gamma; PIASy; Piasg; ZMIZ6; protein inhibitor of activated STAT protein 4; protein inhibitor of activated STAT protein PIASy; protein inhibitor of activated STAT protein gamma; zinc finger, MIZ-type containing 6; protein		

GENE DESCRIPTION	protein inhibitor of activated STAT, 4
GENE ID	51588
ORIGIN SPECIES	Human
ACCESSIONS	BC010047
VECTOR	pDONR223
NCBI SYMBOL	PIAS4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]