



Reference: PTXBC010064

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET59; nicalin; nicalin homolog; nicastrin-like protein		
For research use only.			

GENE DESCRIPTION	nicalin homolog (zebrafish)
GENE ID	56926
ORIGIN SPECIES	Human
ACCESSIONS	BC010064
VECTOR	pDONR223
NCBI SYMBOL	NCLN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]