



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SCPEP1-serine carboxypeptidase 1 Gene

Reference: PTXBC010078

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HSCP1; RISC; retinoid-inducible serine carboxypeptidase; serine carboxypeptidase 1 precursor protein; serine carboxypeptidase 1

For research use only.

GENE DESCRIPTION	serine carboxypeptidase 1
GENE ID	59342
ORIGIN SPECIES	Human
ACCESSIONS	BC010078
VECTOR	pDONR223
NCBI SYMBOL	SCPEP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcctgctcagctcggcgctctcccgctccggctgggtgctcgtcgtcgtccgtcgtcgtcggcgctgaacccagagagctgtcattgactgcgccacagaggaaggcgaaggatgggattatgtagcgtcccgcaaggatgcctacatgtctgtggtcctattatgccacaact
cctgcaagaactctcagaactgcccttggtcattggttcctcaggcggtccaggcggttctagcactggatttggaaacttgggaaattggggccctgacagtgtatcctaaaccaggaaaaccactggttcaggctgccagctcctctatttgggataatccctggggcactgggtt
cagttatgtgaatgtagtgtgctcctatgccaagaactcgtcgtatgggttcctcagacatagttgtctcctgaagacacttctcattgtgcacaaagaattccagacagcttcattctacatttccagagttcctatggaggaataatggcagctgctgcaattgggtctagacgttataaaggccattc
agcgaggggaccatcaagtgcacttgggggggttgccttgggtgattctcgtatcctccctgttgattcgtctcctgggacacttacctgtacagatctgtctctcctgaagacaaatctgacagagctgtctgtagtgaagctgtcagagcaagtactgaatccgtaataaaggcgctcta
cagagagccccacagagctgttggggggaaagcagaagaatcattgaacacagacacagatgggggtgaactctataacatctataaaagcactccacactctacatgaatgggctgagtctatgacacagacacacactagttgtgttttag