

PSMB4-proteasome (prosome, macropain) subunit, beta type, 4 Gene

Reference: PTXBC010098

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HN3; HsN3; PROS-26; PROS26; proteasome subunit beta type-4; 26 kDa prosomal protein; hsBPROS26; macropain beta chain; multicatalytic endopeptidase complex beta chain; proteasome (prosome, macropain) subunit, beta type, 4; proteasome beta chain; proteasome		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, beta type, 4
GENE ID	5692
ORIGIN SPECIES	Human
ACCESSIONS	BC010098
VECTOR	pDONR223
NCBI SYMBOL	PSMB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagcgttttgggtcgcgggtccgacttggcggggggtccggccccaggacagttttaccgattccgctccactcccgattcctcatggtccggcgtctgcactttacagaggtccaatcacgcggaccagaaccccatggtgaccgggacccagtcctcggcggttaagt
cgaggcgaggagtggtgattcccgagacatgctgggacacggcctctggcctgttccgcaacatctctcgattatgcgagtgcaacaacagtaacatgctgggtgcctctggcgactacgctgatttccagtatttgaagcaagttctcgccagatggtgatgaggagctctg
ggagatggacacagctatagctctagagctattcattcgtgctgaccaggccatgtacagccggcgcctcgaagatgaaccccttgggaacacatggtcatcggaggctatgctgagagagagcttccctggttatgtggacatgcttgggtgtagcctatgaagcccttcgctggcc
actggttatggtgcatacttggctcagcctctgctcgagaagttctggagaagcagccagtgctaaagcagaccgagggcccgacttagtagaacgctcatgcgagtgctgtactaccgagatgccgttcttacaaccggttcaaatcgccactgtcaccgaaaaggtgttgaat
agagggaccattgtctacagagaccaatgggatattgcccacatgatcagtggtttgaatga