

DNAJC4-DnaJ (Hsp40) homolog, subfamily C, member 4 Gene

Reference: PTXBC010145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DANJC4; HSPF2; MCG18; dnaJ homolog subfamily C member 4; DnaJ (Hsp40) homolog, subfamily C, member 4; dnaJ-like protein HSPF2; heat shock 40kD protein 2; multiple endocrine neoplasia type 1 candidate protein number 18; DnaJ heat shock protein family (Hsp4		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 4
GENE ID	3338
ORIGIN SPECIES	Human
ACCESSIONS	BC010145
VECTOR	pDONR223
NCBI SYMBOL	DNAJC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgccctactgcccctgcgcctgtgccggctgtggccccgcaacctccctcccggtcctcggagcggccgcccggcagcgggtccagaccagctactattatgaactgttgggggtgcatcctggtgccagcactgaggaagtaaacgagcttcttccaagtcacaagagg
aagggtgaagcagatgcaccttaa