

RPL28-ribosomal protein L28 Gene

Reference: PTXBC010173

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L28; 60S ribosomal protein L28; ribosomal protein L28		

For research use only.

GENE DESCRIPTION	ribosomal protein L28
GENE ID	6158
ORIGIN SPECIES	Human
ACCESSIONS	BC010173
VECTOR	pDONR223
NCBI SYMBOL	RPL28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcgcatctgcaatggatggtcgtcggaactgtccagtttctgatcaagagggaataagcagacctacagcactgagcccaataacttgaaggcccgaattcttccgtacaacggactgattcacgcaagactgtggcgtggagccggcagccgacggcgaaggtgc
gtggtgtcattaaagcgagatccggccagcggagcctgccacctctatgtcggaaccacatcaacaagaatgctcgcgccacgtcagcagcatcagacatgatccgaagaacaagtaccgccccgacctgcgcatggcagccatccgagggccagcgccatctgcg
cagccagaagcctgtgatggtgaagagggaagcggaccgccccaccaagagctcctga