

UBE2A-ubiquitin-conjugating enzyme E2A (RAD6 homolog) Gene

Reference: PTXBC010175

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HHR6A; MRXS30; MRXSN; RAD6A; UBC2; ubiquitin-conjugating enzyme E2 A; E2 ubiquitin-conjugating enzyme A; RAD6 homolog A; ubiquitin carrier protein A; ubiquitin-conjugating enzyme E2A; ubiquitin-protein ligase A; ubiquitin conjugating enzyme E2 A		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2A (RAD6 homolog)
GENE ID	7319
ORIGIN SPECIES	Human
ACCESSIONS	BC010175
VECTOR	pDONR223
NCBI SYMBOL	UBE2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccaccccggtcggcgccctcatcgaggactcaagaggtgcaggagatccagccggagtcagcgggctccgagagaacaacataatggtgtggaacgcgggtcatttcgggctgaaggaccccggttgaggatggaacatttaaaactacaatagaattcact
gaagaatatccaataaaccactacagttgattgtcttaagatgttccatccaaatgtctatgcagatggtatgtatgtctggacatactcagaaccgttggagtcacacatgatgtgtcttcattctaatacctacagtcctgttggatgaaccaatccaatagtcagcaaaa
gccaggctgctcagctgtaccaggagaacaacgggaatatgaaagcgtgttctgcaatagtagaacaagctggcgtgattgttga