

NY-SAR-48-sarcoma antigen NY-SAR-48 Gene

Reference: PTXBC010176

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sarcoma antigen NY-SAR-48; NY-SAR-48; DGT4; HICE1; HAUS augmin-like complex subunit 8; HEC1/NDC80 interacting, centrosome associated 1; HEC1/NDC80-interacting centrosome-associated protein 1; Hec1-interacting and centrosome-associated 1; HAUS augmin like		

For research use only.

GENE DESCRIPTION	sarcoma antigen NY-SAR-48
GENE ID	93323
ORIGIN SPECIES	Human
ACCESSIONS	BC010176
VECTOR	pDONR223
NCBI SYMBOL	NY-SAR-48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagacaactctgctgagttgaaagaaggcagaaaaaatttattaataatgtgaaggagaaggagaagctacagaaaaaggcccacgagctgaagcgcaggcttctctctcagaggaaagcgggagctggcagatgtcctggatcccagatcagatgctcagccccttc
gaggcagtgggccacacgcttcaaggagcaatacaggacatfcggccacggccctggacactaccaggcagcagctgcccgtgaggtccatccactggaggagatgggagcagctcttagacccctgcagcatgaactggtagccactcagcgcctctgggagaacttgatgtt
gggtattcggaagaaatgtgaggtgctgacttactgagcgaactcaaggacgtgacagcgaagaaaggaccttgagctccgaaggagcttggccaggtgctggaactctccgagaggcaagcaagaggcagccttggcaaccaggaagtctgggaagagaccaggggcat
ggcggcccccagccgggtgatttcaatcaagacagtcctgcagagaatctgggggagcaccgaagaacacgcccctgtctgaggagcacaacccgggtgcctcgtcagccccgctcaggccacgttcatagccaagcgaagattttctcaagcagcagcagaagtcctcc
ccctctctctctgcttcaggaggagacttgatga