

YWHAB-tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, beta polypeptide Gene

Reference: PTXBC010352

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GW128; HEL-S-1; HS1; KCIP-1; YWHAA; 14-3-3 protein beta/alpha; 14-3-3 alpha; brain protein 14-3-3, beta isoform; epididymis secretory protein Li 1; protein 1054; protein kinase C inhibitor protein-1; tyrosine 3-monooxygenase/tryptophan 5-monooxygenase act		

For research use only.

GENE DESCRIPTION	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, beta polypeptide
GENE ID	7529
ORIGIN SPECIES	Human
ACCESSIONS	BC010352
VECTOR	pDONR223
NCBI SYMBOL	YWHAB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgacctatctgttccctggataataacctttaagaataatgtctctgagtcaggcgtgtgtgtgctgcatctagtcccaactatttgggaggctgaggcaggagatcgttgagccaggagttaaaagctgcagtgcctgtgtgttgacactgtgaataactgcactccagcctggc
aacatag