

MRPS18A-mitochondrial ribosomal protein S18A Gene

Reference: PTXBC010364

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HumanS18b; MRP-S18-3; MRPS18-3; S18bmt; 28S ribosomal protein S18a, mitochondrial; 28S ribosomal protein S18-3, mitochondrial; MRP-S18-a; S18mt-a; mitochondrial ribosomal protein S18-3; mitochondrial ribosomal protein S18A		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S18A
GENE ID	55168
ORIGIN SPECIES	Human
ACCESSIONS	BC010364
VECTOR	pDONR223
NCBI SYMBOL	MRPS18A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccctcaaggctctggtgctcggtgtggcggtctctcggtggtgctactagcgggcccggcagcgaccagctggtctcggttcagctcggtggtcaggggaagtgggtggagaccaagaagggaagacaactataattgaaggccgtatcacagcgactcccaaggag
agtcctcaatctcttaacccctctggccagtggccatctgcccgttggaaactgaagcacaagtataactatgacgagttctgctgcttagccagttcatccgctcatggagcagctgccccgaagatcacaggcctatgccagggaagaacacgcgaagatcgaggagtgtgtaa
gatggccaccgagcaggtctatfaccaatcacaggcctcggttctgaaggagttgttccgaagagcaaacccaactcaaccggtacctgacgcgtgggtctctggtcctgcaagccatctacaaaaaggccccgctggaacagggtgcgcatgccgtgggtcaccc
cttctgaggacaatgtctgctactcaagaacaccttgaagctgtatcactga