



Reference: PTXBC010370

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	magnesium uptake/transporter TUSC3; D8S1992; M33; MRT22; MRT7; N33; OST3A; tumor suppressor candidate 3; Putative prostate cancer tumor suppressor; oligosaccharyltransferase 3 homolog A		

For research use only.

GENE DESCRIPTION	tumor suppressor candidate 3
GENE ID	7991
ORIGIN SPECIES	Human
ACCESSIONS	BC010370
VECTOR	pDONR223
NCBI SYMBOL	TUSC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]