

TSPAN31-tetraspanin 31 Gene

Reference: PTXBC010377

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tetraspanin-31; sarcoma-amplified sequence; transmembrane 4 protein; tspan-31; tetraspanin 31		

For research use only.

GENE DESCRIPTION	tetraspanin 31
GENE ID	6302
ORIGIN SPECIES	Human
ACCESSIONS	BC010377
VECTOR	pDONR223
NCBI SYMBOL	TSPAN31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtttgtggcctttgctcctcaagaatgcctttgcctcctcaacgtggtctacatgctggtgagctttgtctcattggagtgctgctgggcaaggcctgggtctggtgctcagatccatcatcggcgagtcattgctggtggagctctctctcttattgcagtgctg
gactggtgggtgctgcaaccaccaccaagtcctgctgtcttttacatgatcatccttggtttgctctcctcctcaatttgaatctctgctcctgctgcttataaccgaagcaaacagacagatgtcatcaatgctcttgggtgctgagcaacaagactcgggatgaactgaaaagaa
gttttgattgtggtcttattcaacctcaaccctgatcaacaagattatgattctgactgcaatctgcaagagccagagcccatgccagatgtggagaaaagtttcttaagcattcagacgaagccctgaaaatcttaggggtgtggactctcttagctttacagagatcctg
gtgtttgctagcaatgagatttcggaatcagaagatcctagagccaacccagtcctttctatga