



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TSPAN17-tetraspanin 17 Gene

Reference: PTXBC010405

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBX23; FBXO23; TM4SF17; tetraspanin-17; F-box only protein 23, transmembrane 4 superfamily member 17; tetraspan protein SB134; tetraspanin 17		

For research use only.

GENE DESCRIPTION	tetraspanin 17
GENE ID	26262
ORIGIN SPECIES	Human
ACCESSIONS	BC010405
VECTOR	pDONR223
NCBI SYMBOL	TSPAN17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcgaagcagcacagatcttcaggaacctgaggtcgctgctgcgggaaatctctctgttggcttcacatctgtctcttgggtctggagaccgtcttcctggctatcgcgcctctggcctggggfagagaaggcgctctctgaacatctcagcgtcagacatctgggagcctt
gaccccgctgtggctgtgttgggtagttggaggcgctatcgtggtctcgggcttctcggctcattggggccctccggggagaacacctctctcgaagtcttctcgtgtctcctggtctcatcttctctggagctggcaacggggatctggccttctctcaaggactgattcagacc
agctcaaacctcttcatacaacaacactcgaagcctaccgggagcagatgaccttcaggaacctcattgacttctcaggaatactgcttctcgtcggagcccgaggccccaatgactggaactcaatatctacttcaactgactgacttgaaccccaacggggagcgtcgggggtg
cctcttctcgtcgtcaggggacccctcgggggagctgcctctcaacaccagctgggtcagcagctccgctcaaaactggagcagcagggttcattccacaccaaaaggctcgtctggccagtttgaagaagtgcctcaggacaacctgatttgggtggcgggagttctatggg
catcgccctctccagatcttggcgtcgggagcctccagcctcagcctcagcgaatgaatga