

ENAH-enabled homolog (Drosophila) Gene

Reference: PTXBC010414

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ENA; MENA; NDPP1; protein enabled homolog; mammalian enabled variant 11a; mammalian enabled variant pan; enabled homolog (Drosophila)		

For research use only.

GENE DESCRIPTION	enabled homolog (Drosophila)
GENE ID	55740
ORIGIN SPECIES	Human
ACCESSIONS	BC010414
VECTOR	pDONR223
NCBI SYMBOL	ENAH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaattcaagaagacaactacaagaacagcaacggcaaaaggagctggagcgggaaaggctggagcgagaagaatggaaagagaaggctggagagagaggttagaaagggaaggctggagagggagcgactggacaagaacagctggagagagagagaca
agaacgggaacggcagggaacgcctggagcggcagggaacgcctggagcggcagggaacgcctggatcgggagagggcaagaagaacagagagagcctggagagactggaaacgggagagggcaagaaggagcgacaagagcagttagaagggaacagctggaatg
ggagagagagcgagaaatacaagtgcgtgtaagaagttaacaactctctgtgtcttctgtatctgagcttgatgtctctgtgcagtcattgttaa