



Reference: PTXBC010446

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf4; transmembrane protein 59-like; brain-specific membrane-anchored protein; transmembrane protein 59 like		

For research use only.

GENE DESCRIPTION	transmembrane protein 59-like
GENE ID	25789
ORIGIN SPECIES	Human
ACCESSIONS	BC010446
VECTOR	pDONR223
NCBI SYMBOL	TMEM59L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]