



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD14-CD14 molecule Gene

Reference: PTXBC010507

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD14 molecule; monocyte differentiation antigen CD14; myeloid cell-specific leucine-rich glycoprotein		

For research use only.

GENE DESCRIPTION	CD14 molecule
GENE ID	929
ORIGIN SPECIES	Human
ACCESSIONS	BC010507
VECTOR	pDONR223
NCBI SYMBOL	CD14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr