



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FLJ22167-hypothetical protein FLJ22167 Gene

Reference: PTXBC010609

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	I-GlcNAc-6-ST; I-GlcNAc6ST; glcNAc6ST-3; gn6st-3; hIGN6ST; carbohydrate sulfotransferase 5; GST4-alpha; N-acetylglucosamine 6-O-sulfotransferase 3; carbohydrate (N-acetylglucosamine 6-O) sulfotransferase 5; galactose/N-acetylglucosamine/N-acetylglucosamin		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ22167
GENE ID	79583
ORIGIN SPECIES	Human
ACCESSIONS	BC010609
VECTOR	pDONR223
NCBI SYMBOL	FLJ22167
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr