



Reference: PTXBC010623

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LCR-F1; NRF1; TCF11; nuclear factor erythroid 2-related factor 1; NF-E2-related factor 1; NFE2-related factor 1; TCF-11; locus control region-factor 1; nuclear factor, erythroid derived 2, like 1; transcription factor 11 (basic leucine zipper type); trans		

For research use only.

GENE DESCRIPTION	nuclear factor (erythroid-derived 2)-like 1
GENE ID	4779
ORIGIN SPECIES	Human
ACCESSIONS	BC010623
VECTOR	pDONR223
NCBI SYMBOL	NFE2L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]