

ZNF22-zinc finger protein 22 (KOX 15) Gene

Reference: PTXBC010642

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HKR-T1; KOX15; ZNF422; Zfp422; zinc finger protein 22; krox-26 protein; zinc finger protein 22 (KOX 15); zinc finger protein KOX15; zinc finger protein Krox-26		

For research use only.

GENE DESCRIPTION	zinc finger protein 22 (KOX 15)
GENE ID	7570
ORIGIN SPECIES	Human
ACCESSIONS	BC010642
VECTOR	pDONR223
NCBI SYMBOL	ZNF22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggttagcaaaagcctaaagcgggtattctcggagctcaagcgaagcctatgagaacaagcgaagcagccggcagcggcagaagtgggcgatgactattcgattgactcaagctcagtagactcagaagaagcttgatgacaaacctataaatgtactgaatgt
gaaaagagtttcagtcagagttcaactcttttcaaccagagaagatccataactggaagaaatccataaatgtgctgattgtgggaaagtgtcttcagagtttaattcattcagcatcgacggatccatagggggaagcctacaaatgtgatgagtgaggagaagcttcaaacag
agctcaaatcattcagcaccagagaattcactggagaaaaacctatcagtgatgagtggtggcgggtgttcagccagagctccaccctattcaacatcagagaaccacactggggagaaccctaccagtgcatgtgaatgtggcaaatgttcagtcagagctctcatctgagg
cagcacatgaaggtgcataaagaagaagcctgtaaacccggggcaaaatcagggtgaagactcactaccctctggaaagctgtacaggaaggaagctgtggctgtctccgttaa