



Reference: PTXBC010738

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf26; HMFN1876; JFP11; TMEM112; TMEM112A; lipase maturation factor 1; transmembrane protein 112		

For research use only.

GENE DESCRIPTION	lipase maturation factor 1
GENE ID	64788
ORIGIN SPECIES	Human
ACCESSIONS	BC010738
VECTOR	pDONR223
NCBI SYMBOL	LMF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]