

DYNLL2-dynein, light chain, LC8-type 2 Gene

Reference: PTXBC010744

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNCL1B; Dlc2; RSPH22; dynein light chain 2, cytoplasmic; 8 kDa dynein light chain b; DLC8b; radial spoke 22 homolog; dynein light chain LC8-type 2		

For research use only.

GENE DESCRIPTION	dynein, light chain, LC8-type 2
GENE ID	140735
ORIGIN SPECIES	Human
ACCESSIONS	BC010744
VECTOR	pDONR223
NCBI SYMBOL	DYNLL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgaccggaaggcagtgatcaagaacgcagacatgtctgaggacatgcaacaggatgccgttgactgcgccacgcaggccatggagaagtacaatatagagaaggacattgctgcctatatcaagaaggaattgacaagaaatataacacctacctggcattgtatcgtggccga
aatgttgagctacgtcacacacagacaaagcacttcattctatcttttacttgggtcaagttgcaatcctcctcttcaagtcaggctag