

SNX24-sorting nexin 24 Gene

Reference: PTXBC010886

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO1284; SBBI31; sorting nexin-24; sorting nexin 24; sorting nexin 24		
For research use only.			

GENE DESCRIPTION	sorting nexin 24
GENE ID	28966
ORIGIN SPECIES	Human
ACCESSIONS	BC010886
VECTOR	pDONR223
NCBI SYMBOL	SNX24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtctacatccgctcttgcctatgaagagagcgacctggagcggggatcacggtgtttaagatagaagtgcctaatgaatggagaaaacattttgtgaaagagatacagcgaatttcagctttgcacaaaaagcctaagaaatgtataaaactccagaatcccttctaaaca
tgftaggaaactgggtcccaaaagtcttgaacagcgcgacgacaaggcttggaaacatactacaggctgtcatttagaaaatgaagaacttcccaactgtttcttgatttcctaaatgtgcgacacttgcctctctaccaaaaggcagaagtgtggatctttgatgaacagagctgaagag
tcaagcaaacgtgccaccagcctgtgctgtgtcctcaggatccatgtcttgcctgcagcagcgattttccaaatgtgtgtattgaaggagtcctccatgggatattttaccctcatctacagcccaggtag