

IER3IP1-immediate early response 3 interacting protein 1 Gene

Reference: PTXBC010888

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC039; MEDS; PRO2309; immediate early response 3-interacting protein 1; immediate early response 3 interacting protein 1		

For research use only.

GENE DESCRIPTION	immediate early response 3 interacting protein 1
GENE ID	51124
ORIGIN SPECIES	Human
ACCESSIONS	BC010888
VECTOR	pDONR223
NCBI SYMBOL	IER3IP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctttaccctgtactcactgctgcaggcagccctgctctgcgtcaacgccatcgagtgctgcacgaggagcgttctcaagaacattgctggggaacagaccagggaattgggtgattggagaagagccgggaattaaatcacagctaatgaacctattcgatctgtaagaac
cgtgatgagagtgccattgataatagtaaacatgcaattgtgtactttattatttgatga