

DERL2-Der1-like domain family, member 2 Gene

Reference: PTXBC010890

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-101; DERtrin-2; F-LAN-1; F-LANa; FLANa; derlin-2; Der1-like domain family, member 2; carcinoma related; degradation in endoplasmic reticulum protein 2; derlin 2		

For research use only.

GENE DESCRIPTION	Der1-like domain family, member 2
GENE ID	51009
ORIGIN SPECIES	Human
ACCESSIONS	BC010890
VECTOR	pDONR223
NCBI SYMBOL	DERL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtaccagagcttgcggctggagtagctgcagatccaccggcagcggcctacaccactgcctgcgtcctaccaccggcggcgtgcagttggaattgacacacctttcagttgtacttcaatcctgaattaatctttaaactttcaaatatggagattaatcaccactcttatt
ttttggccagttggattcaattttttatatacatgattttctatctggtactgtcgaatgctagaagaaggctcttccgaggtcggacagcagactttgtattatgttcccttttgggtgattctaatgacctttttgggtctgtttgtgagcttagttttctggccagggcctttacaataatgctcgtct
atgtgtggagccgaaggaaacctatgtccgatgaacttctcggccttcaacttccaggcccccttctgcccctgggtgctcatgggaatttctctgtttgtgggaactcaatcattgtggaccttttgggtattgcagttggacacatatattttcttggagatgtatttccaatcaacctg
gtggaataagaattctgaaaacacatctattttgaaagctattttgataccagatgaggatccaaattacaatccactacctgaggaacggcaggaggctcgcctgggtgagggcagcggcttgagggttaa