



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CHMP4A-chromatin modifying protein 4A Gene

Reference: PTXBC010893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf123; CHMP4; CHMP4B; HSPC134; SHAX2; SNF7; SNF7-1; VPS32-1; VPS32A; charged multivesicular body protein 4a; SNF7 homolog associated with Alix-2; Snf7 homologue associated with Alix 2; chromatin modifying protein 4A; vacuolar protein sorting-associate		

For research use only.

GENE DESCRIPTION	chromatin modifying protein 4A
GENE ID	29082
ORIGIN SPECIES	Human
ACCESSIONS	BC010893
VECTOR	pDONR223
NCBI SYMBOL	CHMP4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr