

PAGE4-P antigen family, member 4 (prostate associated) Gene

Reference: PTXBC010897

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT16.7; GAGE-9; GAGEC1; JM-27; JM27; PAGE-1; PAGE-4; P antigen family member 4; P antigen family, member 4 (prostate associated); g antigen family C member 1; prostate-associated gene protein 4; PAGE family member 4		

For research use only.

GENE DESCRIPTION	P antigen family, member 4 (prostate associated)
GENE ID	9506
ORIGIN SPECIES	Human
ACCESSIONS	BC010897
VECTOR	pDONR223
NCBI SYMBOL	PAGE4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcacgagtgagatcaagatccagaggagaggagatgggtcaggaggctcccgatgtgttcattcgtgctcccggtgaatctcagcaagaggaaccacaaactgacaatcaggatattgaacctggacaagagagagaaggaaacacctccgacgaagaacgtaaagta
gaagggtattgccaggaaatggatctgaaaagactcggagtgagcgtggagatggctctgatgtaaaagagaagactccacctaatacctaagcatgctaagactaaagaagcaggagatgggcagccataa