

C6orf203-chromosome 6 open reading frame 203 Gene

Reference: PTXBC010899

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf203; HSPC230; PRED31; chromosome 6 open reading frame 203		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 203
GENE ID	51250
ORIGIN SPECIES	Human
ACCESSIONS	BC010899
VECTOR	pDONR223
NCBI SYMBOL	C6orf203
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctatggctagtgtaaatctgctccgggtgtttaaagaagccagatgcctggattggactctgggggtgttctccgaggacaccttcatacaaaactctgtactctctggaatcgatactgtattttctagtagcaagttacgtgcaccaattataaaacactttttataatattttctact
gagactcccagggtctttactatctccagaatgtattttctttccgtaagactcaaaagtataataaggtctacaaaactctactaaaaagtctctgcaaaaagtagatgaagaggactctgatgaagaaagccatcatgatgagatgagtgagcagggaagaggagcttgaggatgacctac
ttagtcaaaaactataaagacctggaaaaagcagttcagcttttcgggtatgatgttctctgaagacggggctagatattgggagaacaagtggaagatgctttctacaaaggtgaactcaggctgaatgaggaaaaattatggagaaaaagcagaacgggtgaaagtgggagatacatt
ggatctctcattggagagataaagaagcagggaacagacagattatgcggattctctgaaaaaagtgtttgaagagaagactgaaagtgaataacagagtggtgttacggcggtggaaggtttaagttgcctaagaagagaatgtctaataa